CozyLine / CozyTube

- DE** Beheizter Atemschlauch
- EN** Heated Breathing tube
- FR** Circuit respiratoire chauffant
- IT** Tubo di respirazione riscaldato
- NL** Verwarmde beademingslang

Zeichenerklärung
Key to symbols
 Signification des symboles
Descrizione dei simboli
 Uitleg van symbolen

Gebrauchsanweisung befolgen
Follow the instructions
 Se reporter au mode d'emploi
Rispettare le istruzioni d'uso
 Volg de gebruiksaanwijzing

Eindeutige Produktkennzeichnung
Unique Device Identifier
 Identificateur de produit unique
Identificazione unica del dispositivo
 Unieke apparaatidentificatie

Chargennummer
Lot number
 Numéro de lot
Numero di lotto
 Batchnummer

Hersteller
Manufacturer
 Producteur
Produttore
 Fabrikant

Hersteldatum
Date of manufacture
 Date de fabrication
Data di fabbricazione
 Productiedatum

Verfallsdatum
Use by date
 À utiliser avant
Data di scadenza
 Uiterste gebruiksdatum

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Intended use

WILAMED heated breathing tubes are intended for channelling humidified breathing gas and are an accessory for PAP devices as well as non-invasive patient interfaces with standardized connections.

Medical indication and contraindications

Heated breathing tubes are indicated for patients requiring permanently a positive airway pressure ventilation or mask ventilation. They are used as part of non-invasive respiratory therapy in both the clinical and non-clinical setting. There are no known contraindications.

Intended patient group

The product can be used in children and adults over 30 kg body weight.

User qualifications

The person operating the device is referred to as the user in these instructions for use. Patients can also be users if they meet the necessary prerequisites for operating the product. The device may only be used on the instructions of a doctor.

As an operator or user, you must be familiar with the operation of this medical device. The operator is responsible for ensuring the compatibility of the device and all components or accessories connected to the patient before use.

When handing over the device to the patient, you, as the attending physician or healthcare professional, must carry out a briefing on the function of the device.

Reporting Obligation

All serious incidents that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Disposal

To preserve and protect the environment, to prevent environmental pollution, and to enable recycling of raw materials, the European Commission has determined that electrical and electronic equipment shall be taken back by the manufacturer in order to enable their proper disposal. Devices with the symbol "Not for disposal as household waste" may not be disposed in the unsorted municipal waste. Please note that both the breathing tube and the adapter must be reprocessed in accordance with the descriptions in these instructions for use before they are disposed of.

CozyLine – Use guide

1. Remove the heated breathing tube and the AC adapter from the packaging.
2. Check the product for damage and contamination. Damaged and/or contaminated products may not be used! For the WILAMED product "Heated breathing tube", we recommend cleaning the tube before first use.

WARNING: If the product is combined with other products, the instructions for their use from the respective manufacturers must be respected.

3. Connect the end of the breathing tube where the connection cable is located to the output of the PAP device.
4. Connect the mask-side end of the breathing tube to the mask.
5. Connect the electrical connection cable of the breathing tube to the AC adapter included.

WARNING:

It is important to ensure that the adapter is positioned next to the PAP device in a fall-proof manner. It is not permitted to leave the adapter hanging freely on the connection cable of the breathing tube.

6. Connect the adapter to the power cable supplied and plug it into the mains socket.
7. Disconnect the mains plug from the socket after use.

CozyTube – Use guide

1. Remove the heated breathing tube from the packaging.
2. Check the product for damage and contamination. Damaged and/or contaminated products must not be used! For the WILAMED product "Heated breathing tube", we recommend cleaning the tube before first use.

WARNING: If the product is combined with other products, the instructions for their use from the respective manufacturers must be respected.

3. Connect the end of the breathing tube where the connection cable is located to the output of the PAP device.
4. Connect the mask-side end of the breathing tube to the mask.
5. Connect the electrical connection cable of the breathing tube to the connection socket of the CozyLine adapter.

Warnings

Failure to observe the following warnings may affect product performance or safety:

- (1) Store product in a cool, dry place away from direct sunlight and do not store or operate near heat sources.
- (2) Before use, make sure the heated breathing tube is suitable for the intended respiratory therapy device.
- (3) Observe operating instructions for the respiratory therapy device and accessories.
- (4) This product may only be connected by persons who have been instructed by trained specialist personnel.
- (5) This product is not suitable for the supply of flammable gases or medication and not for inhalation purposes.
- (6) Risk of burns:
 - a Do not cover heated breathing tubes with insulating material (e.g. cloths, blankets, etc.).
 - b To avoid overheating, a continuous, minimal flow of breathing gas in the heated breathing tube is required.
 - c Patient interface must not be locked or blocked.
- (7) Heated breathing tubes are only to be used in conjunction with compatible therapy devices; otherwise the quality of the therapy may be impaired.
- (8) Do not use the product at room temperatures above 30°C, to avoid elevated temperatures of the air inhaled.
- (9) To connect and disconnect the product from the therapy device and/or mask, always grip the grip surfaces of the tube sleeve to prevent damage to the breathing tube.
- (10) The product must not be used in rooms where EM disturbances of high intensity occur (e.g. MRI). In addition, the function of the product may be impaired by RF surgical equipment. Otherwise, the performance of this device may be reduced.
- (11) Never place the product around the neck, to avoid strangulation.
- (12) Danger of suffocation due to the use of mouth-nose masks without exhalation valve! Please read the user manual of the mouth-nose masks before starting therapy.
- (13) Prevent the tube or tube system from loosening during use, especially during outpatient use.
- (14) Do not bend the tube or cable.
- (15) Do not use the product if the insulation or other parts are damaged.
- (16) Do not connect the product to incompatible other products or accessories.
- (17) Maintenance must not be carried out while the device is in use.
- (18) This device may not be changed without the manufacturer's permission, otherwise the quality of the therapy may be impaired or the patient may be injured.
- (19) Keep therapy equipment and accessories away from children and pets.
- (20) For the compatibility of this product with other products, the following must be observed: WILAmEd's heated breathing tubes have been developed in accordance with the applicable standards EN ISO 5367 and EN ISO 5356-1. Therefore, heated breathing tubes may be used with other ventilators or respiratory therapy devices, provided that the following requirements are met:
 - a The connection points of the corresponding PAP respiratory therapy devices are standardized according to applicable standards EN ISO 5367 and EN ISO 5356-1.
 - b Before connecting the patient, always perform all pre-commissioning checks, such as a tube system test, to avoid possible injury to the patient. Ensure that an alternative therapy option is available.
 - c The heated breathing tubes meet the requirements for the corresponding PAP respiratory therapy devices with regard to general safety and performance.
 - d There are no known contradictions with regard to the intended use, the intended patient population, the intended user or the intended application environment.
 - e CozyLine: The responsible organization has taken into account that by combining the product with a respiratory therapy device, an ME system is established. In doing so, the requirements for ME systems from the latest edition of IEC 60601-1 must be observed.
- (21) CozyLine: The use of accessories, transducers and cables other than those specified herein or supplied by the manufacturer of this equipment may result in increased electromagnetic emissions or reduced electromagnetic immunity to interference of this equipment and may result in malfunctions.
- (22) CozyLine: Portable RF communication equipment (such as mobile phones, including peripheral parts such as antenna cables and external antennas) should not be used at a distance of less than 30 cm from any part of the product, including the cables specified by the manufacturer. Otherwise, the performance of this device may be reduced.
- (23) CozyLine: The product complies with the applicable limits for radiated RF emissions to avoid adverse events in the patient and user due to electromagnetic interference. In order to avoid adverse events in the patient and user due to electromagnetic interference, the product must not be modified or altered.
- (24) CozyLine: Avoid using this device in the immediate vicinity of or stacked with other equipment, as this can result in faulty operation. However, if it is necessary to use this equipment in the manner described above, this and other devices must be monitored to ensure their proper operation.
- (25) CozyLine: For fire safety and other reasons, the use of multiple sockets is not recommended. The responsible organization is advised that connecting electrical equipment to a multiple socket will result in the creation of an ME system and may lead to a reduced level of safety.
- (26) CozyLine: Only use the AC adapter supplied. Refer to the specifications of the adapter in the "Technical data" section.
- (27) CozyLine: To minimize the risk of overheating, unplug the power cable after use.

Reprocessing

WARNINGS:

- (1) Reusable products must be reprocessed after use, as otherwise there is an increased risk of infection. Weekly reprocessing is recommended. The local hygiene and reprocessing regulations take priority and must be followed.
- (2) The manufacturer's instructions for the cleaning agents, disinfectants and reprocessing devices must be observed, as well as the manufacturer's instructions for the reprocessing of the accessories.
- (3) Do not use ozone cleaning equipment.
- (4) The electrical connections must not be immersed in liquid and no liquid should get into the connections, as otherwise there is a risk of electric shock.
- (5) With reprocessed products, signs of wear may occur, e.g. cracks, deformations, discoloration or detachment. Check products for signs of wear and replace if necessary.

Reprocessing cycles

Breathing tube:

Weekly reprocessing is recommended. Reprocessing the heated breathing tube is limited to a maximum of 30 reprocessing cycles. If reprocessing is no longer possible or if signs of wear occur, WILAmEd recommends replacing the product.

AC adapter:

The cleaning of the adapter can be carried out as required within the specified service life. To avoid power reduction or shutdown due to overheating, dust and/or dirt deposits must be removed regularly.

Grading and classification

The following classification is according to Spaulding, and is a WILAmEd recommendation:

Part	Classification	Reprocessing
Breathing tube heated	non-critical	Cleaning, disinfection

Validated reprocessing instruction

Initial processing at the place of use	Initial preparation at the place of use is recommended.	
Preparation before cleaning	Breathing tube	Disassemble the heated breathing tube and allow to cool completely.
	AC adapter	The adapter must be disconnected from the mains by unplugging it. The breathing tube must be disconnected.
Manual cleaning	Breathing tube	<ol style="list-style-type: none"> 1. Preparation of cleaning immersion bath: Detergent: - WILAsil or Löwosan: 2 teaspoons concentrate (approx. 10 ml)/2 to 3 l, 25 °C warm tap water 2. Soak the breathing tube in cleaning solution for 3 minutes! <p>WARNING: Electrical connections must remain dry!</p> <ol style="list-style-type: none"> 3. Ensure all surfaces are fully moistened. 4. Brush the breathing tube with a soft brush for at least 30 seconds to remove coarse contaminants. 5. Carefully remove the heated breathing tube from the cleaning solution. 6. Rinse the heated breathing tube under clean, running water until no residues are visible. To do so, hold one end of the breathing tube closed and fill the other with clean, running water. Then also close the other end of the tube and move it vigorously back and forth. Drain the water from the breathing tube again. Repeat at least 5 times until no more detergent residues are visible. 7. Allow the breathing tube to dry well, for example by hanging it over a clothesline. <p>WARNING: Avoid direct sunlight!</p>
	AC adapter	<ol style="list-style-type: none"> 1. Use a soft, lint-free cloth to wipe down the power supply housing and the power supply and tube connection cable, including the plugs. <p>WARNING: Electrical connections/plugs must remain dry! No moisture should be allowed to penetrate the openings.</p> <p>WARNING: Avoid direct sunlight!</p>

Manual disinfection	Breathing tube	<ol style="list-style-type: none"> 1. Actively fill the disinfectant ASP CIDEX OPA Solution #20391 into the disinfectant tank (Cidex OPA). WARNING. The solution is a ready-to-use solution and the service life must be checked using the CIDEX OPA solution test strips before each use. 2. Allow the breathing tube to soak completely in the disinfectant solution for 5 minutes. WARNING: Electrical connections must remain dry! 3. Ensure all surfaces of the breathing tube are fully moistened. 4. After the exposure time, remove the breathing tube from the disinfectant. 5. To rinse: Allow the breathing tube to soak completely in a tank of deionized water for at least 1 minute. WARNING: Electrical connections must remain dry! 6. Repeat Step 5 twice to remove the disinfectant solution completely. 7. Remove the breathing tube from the tank and let it dry well, for example by hanging it over a clothesline. WARNING: Avoid direct sunlight!
Maintenance, control and inspection	Breathing tube	Check the heated breathing tube for visible damage such as cracking, softening or hardening, and replace the product if necessary.
	AC adapter	Check the adapter, power supply and tube connection cables for visible damage such as damaged adapter housing, cracking, chafing and exposed wires and replace if necessary.
Packaging	The products may be stored in their original packaging until reuse.	
Storage	The products must be stored dry and free of contamination in the original packaging. Re-contamination and damage during transport must be avoided. The information accompanying the product regarding storage and transport must be observed.	

The instructions listed above have been validated by the medical device manufacturer as suitable for preparing the medical device for reuse. The person reprocessing the device is responsible for ensuring that the actual reprocessing carried out with the equipment, materials and personnel used in the reprocessing facility achieves the desired result. This requires verification and/or validation and routine monitoring of the procedure.

Other material-compatible resources

- pH-neutral detergent
- VAH-listed high-level disinfectant solution based on ortho-phthalaldehyde

Accessories and options

WILAmEd has neither approved accessories nor options for the heated breathing tube product.

Maintenance and spare parts

When used as intended, the breathing tube as well as the adapter and power cable are maintenance-free within the specified service life. The user must check the breathing tube as well as the adapter and power cable for visible damage such as cracking, softening or hardening before use. The adapter must not have any damage to the housing, tube connection cable, or power cable. If necessary, the respective product must be replaced.

Operators must comply with local regulations on periodic inspections. It is recommended that an electrical safety test be carried out annually in accordance with IEC 62353.

The breathing tube should not be used for more than 12 months and must be replaced after this time. The included adapter has a 5-year service life. It must be replaced after that time. Contact your doctor, respiratory therapist or nurse for assistance.

CozyLine - Troubleshooting

Fault	Cause	Elimination
Breathing tube does not warm up	No power supply available	Check the adapter power cable for a secure connection Make sure the plug of the breathing tube is undamaged and connected to the adapter
	Adapter defective	Contact your doctor, respiratory therapist, nurse, or caregiver.
	Heating wire in the breathing tube defective	Replace breathing tube

CozyTube - Troubleshooting

Fault	Cause	Elimination
Breathing tube does not warm up	No power supply available	Make sure the plug of the breathing tube is undamaged and connected to the adapter.
	Heating wire in the breathing tube defective	Replace breathing tube

Technical specifications

CozyLine – Product specifications

Supply voltage	100 – 240 VAC
Mains frequency	50/60Hz
Power consumption	max. 40 VA
Protection class	II
Maximum temperature of the released gas	< 43°C

CozyTube – Product specifications

Supply voltage		12,0 – 12,5 V DC
Electrical resistance	2700621 - BT2634 (ø15mm)	12,8 ± 1,0 Ω
	2700631 - BT2630 (ø19mm)	11,4 ± 1,1 Ω

Technical description

The heated breathing tube with external AC adapter can be used with all therapy devices for which there is a declaration of compatibility. The validation was conducted with the prisma SMART max therapy device (WM090TD). A list of all compatible therapy devices is available on request.

The heated breathing tube is an active medical device that receives power for its heating wire from the adapter included. The heated breathing tube is considered a type BF applied part, providing additional protection for the patient against electrical hazards. The heated breathing tube has no essential performance features or measuring function, but transports the heated and humidified breathing gas from the therapy device to the patient.

Operating conditions

Ambient temperature: +5 °C to +30 °C
Relative humidity: 10 % to 90 %, non-condensing
Ambient pressure: 700 hPa to 1060 hPa

Environmental conditions during transport and storage

Ambient temperature: -20 °C to +60 °C
Relative humidity: 10 % to 90 %, non-condensing
Ambient pressure: 700 hPa to 1060 hPa

CozyTube & CozyLine – Breathing tube information

Supply voltage	12.0 – 12.5 VDC
Maximum limit pressure	60 hPa
Maximum flow rate, 15mm	150 l/min
Maximum flow rate, 19 mm	170 l/min
Inspiratory resistance, 15mm	≤ 0,12 hPa / (l/min) @ 15 l/min
Compliance, 15mm	≤ 4 ml/hPa @ 60±3 hPa
Inspiratory resistance, 19mm	≤ 0,06 hPa / (l/min) @ 30 l/min
Compliance, 19mm	≤ 5 ml/hPa @ 60±3 hPa

CozyLine – Adapter information

Manufacturer	DONGGUAN SHILONG FUHUA ELECTRONIC CO LTD
Model No.	UES24D1-120200SPA / UE220914STEJ1RM
Supply voltage	100 – 240 VAC
Mains frequency	50/60Hz
Output voltage	12.0 V

CozyLine – Electromagnetic interference

Interference test	Conformity
Conducted and radiated emissions (CISPR 11)	Group 1 / Class B
Distortion due to harmonic currents (IEC 61000-3-2)	Class A
Voltage fluctuations and flicker (IEC 61000-3-3)	compliant

CozyLine – Electromagnetic immunity to interference

Interference immunity testing	Interference immunity test levels
Electrostatic discharge (IEC 61000-4-2)	± 8 kV contact / ± 15 kV air
Radio-frequency electromagnetic fields (IEC 61000-4-3)	10V/m (80 MHz – 2.7 GHz)
Radio-frequency electromagnetic fields from wireless communication devices (IEC 61000-4-3)	9 – 28 V/m (385 MHz – 5.785 GHz)
Electrical fast transient disturbances/bursts (IEC 61000-4-4)	± 2 kV for grid connection
Surge voltages (IEC 61000-4-5)	± 1 kV line to line
Conducted disturbances induced by radio-frequency fields (IEC 61000-4-6)	3 Vrms (150 kHz – 80 MHz) 6 Vrms (in ISM and amateur radio frequency bands between 150 kHz and 80 MHz)
Power frequency magnetic fields (50/60 Hz) (IEC 61000-4-8)	30A/m
Voltage dips and voltage interruptions (IEC 61000-4-11)	0% UT; 1/2 period / 0% UT; 1 period / 70% UT; 25/30 periods 0% UT; 250/300 periods
Close proximity magnetic fields (IEC 61000-4-39)	8 A/m @ 30 kHz / 65 A/m @ 134.2 kHz / 7.5 A/m @ 13.56 MHz

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